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



BIOFIRE® FILMARRAY® TORCH*

**BIOFIRE® FILMARRAY®
Panel**

* Calculations based on running the BIOFIRE® FILMARRAY® Respiratory 2 (RP2) Panel over a 24 hour day.

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
**Up to 22
Targets**

- 3-4 bacteria
- 11–18 viruses

**BIOFIRE®
BLOOD CULTURE
IDENTIFICATION 2
PANEL†***

**43
Targets**

- 26 bacteria
- 7 yeast
- 10 antimicrobial resistance genes

**BIOFIRE®
FILMARRAY®
GASTROINTESTINAL
PANEL†‡**

**22
Targets**


- 13 bacteria
- 5 viruses
- 4 parasites

**BIOFIRE®
FILMARRAY®
MENINGITIS/
ENCEPHALITIS
PANEL†***

**14
Targets**

- 6 bacteria
- 7 viruses
- 1 yeast

**BIOFIRE®
FILMARRAY®
PNEUMONIA
PANEL†***

**33
Targets**

- 18 bacteria
- 8 viruses
- 7 antimicrobial resistance genes

**BIOFIRE®
JOINT
INFECTION
PANEL†***

**39
Targets**

- 29 bacteria
- 2 yeast
- 8 antimicrobial resistance genes

Pathogen target number varies between the BIOFIRE RP2.1†, RP-EZ, and RP2.1-EZ (EUA)* Panels.

BIOFIRE® RESPIRATORY 2.1 PLUS (RP2.1plus) PANEL

23 Targets. ~45 Minutes.

VIRUSES

Adenovirus
Coronavirus 229E
Coronavirus HKU1
Coronavirus NL63
Coronavirus OC43
Middle East respiratory syndrome coronavirus (MERS-CoV)
Severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2)
Human metapneumovirus
Human rhinovirus/enterovirus

Influenza A virus
Influenza A virus A/H1
Influenza A virus A/H3
Influenza A virus A/H1-2009
Influenza B virus
Parainfluenza virus 1
Parainfluenza virus 2
Parainfluenza virus 3
Parainfluenza virus 4
Respiratory syncytial virus

BACTERIA

Bordetella parapertussis
Bordetella pertussis
Chlamydia pneumoniae
Mycoplasma pneumoniae

BIOFIRE RP2.1plus Panel Specifications

Sample Type: nasopharyngeal swab in transport media or saline

Sample Volume: 0.3 mL

Storage Conditions: all kit components stored at room temperature (15-25 °C)

Hands-On Time: approximately 2 minutes

Overall Performance: 97.1% sensitivity, 99.3% specificity¹

SARS-CoV-2 Performance: 98.4% PPA, 98.9% NPA²

Part Number

BIOFIRE RP2.1plus Panel Reagent Kit (30 Pouches): 423740

BIOFIRE® FILMARRAY® PNEUMONIA *PLUS* (PN*plus*) PANEL

34 Targets. ~1 Hour.

BACTERIA (Semi-Quantitative)

Acinetobacter calcoaceticus-baumannii complex
Enterobacter cloacae complex
Escherichia coli
Haemophilus influenzae
Klebsiella aerogenes
Klebsiella oxytoca
Klebsiella pneumoniae group
Moraxella catarrhalis
Proteus spp.
Pseudomonas aeruginosa
Serratia marcescens
Staphylococcus aureus
Streptococcus agalactiae
Streptococcus pneumoniae
Streptococcus pyogenes

ATYPICAL BACTERIA (Qualitative)

Chlamydia pneumoniae
Legionella pneumophila
Mycoplasma pneumoniae

VIRUSES

Adenovirus
Coronavirus
Human metapneumovirus
Human rhinovirus/enterovirus
Influenza A virus
Influenza B virus
Middle East respiratory syndrome coronavirus (MERS-CoV)
Parainfluenza virus
Respiratory syncytial virus

ANTIMICROBIAL RESISTANCE GENES

Carbapenemases
IMP
KPC
NDM
OXA-48-like
VIM

ESBL
CTX-M

Methicillin Resistance
mecA/C and MREJ (MRSA)

BIOFIRE PN*plus* Panel Specifications

Sample Type: BAL-like (including mini-BAL), Sputum-like (including ETA)

Hands-On Time: approximately 2 minutes

Performance: BAL-like (including mini-BAL): 96.2% sensitivity, 98.4% specificity; Sputum-like (including ETA): 96.3% sensitivity, 97.3% specificity⁶

Storage Conditions: all kit components stored at room temperature (15-25 °C)

Part Number

BIOFIRE PN*plus* Panel Reagent Kit (30 Pouches): RFIT-ASY-0143

BIOFIRE® BLOOD CULTURE IDENTIFICATION 2 (BCID2) PANEL

43 Targets. ~1 Hour.

GRAM-NEGATIVE BACTERIA

*Acinetobacter calcoaceticus-
baumannii* complex
Bacteroides fragilis
Enterobacteriales
Enterobacter cloacae complex
Escherichia coli
Klebsiella aerogenes
Klebsiella oxytoca
Klebsiella pneumoniae group
Proteus spp.
Salmonella spp.
Serratia marcescens
Haemophilus influenzae
Neisseria meningitidis
Pseudomonas aeruginosa
Stenotrophomonas maltophilia

GRAM-POSITIVE BACTERIA

Enterococcus faecalis
Enterococcus faecium
Listeria monocytogenes
Staphylococcus spp.
Staphylococcus aureus
Staphylococcus epidermidis
Staphylococcus lugdunensis
Streptococcus spp.
Streptococcus agalactiae
Streptococcus pneumoniae
Streptococcus pyogenes

YEAST

Candida albicans
Candida auris
Candida glabrata
Candida krusei
Candida parapsilosis
Candida tropicalis
Cryptococcus (C. neoformans/C. gattii)

ANTIMICROBIAL RESISTANCE GENES

Carbapenemases

IMP
KPC
OXA-48-like
NDM
VIM

Colistin Resistance

mcr-1

ESBL

CTX-M

Methicillin Resistance

mecA/C
mecA/C and MREJ (MRSA)

Vancomycin Resistance

vanA/B

BIOFIRE BCID2 Panel Specifications

Sample Type: positive blood culture

Sample Volume: 0.2 mL

Storage Conditions: all kit components stored at room temperature (15-25 °C)

Hands-On Time: approximately 2 minutes

Performance: 99.0% sensitivity, 99.8% specificity³

Part Number

BIOFIRE BCID2 Panel Reagent Kit (30 Pouches): RFIT-ASY-0147

US FDA-cleared |  2797

BIOFIRE® FILMARRAY® MENINGITIS/ENCEPHALITIS (ME) PANEL

14 Targets. ~1 Hour.

BACTERIA

Escherichia coli K1
Haemophilus influenzae
Listeria monocytogenes
Neisseria meningitidis
Streptococcus agalactiae
Streptococcus pneumoniae

VIRUSES

Cytomegalovirus (CMV)
Enterovirus (EV)
Herpes simplex virus 1 (HSV-1)
Herpes simplex virus 2 (HSV-2)
Human herpesvirus 6 (HHV-6)
Human parechovirus (HPeV)
Varicella zoster virus (VZV)

YEAST

Cryptococcus (C. neoformans/C. gattii)

BIOFIRE ME Panel Specifications

Sample Type: cerebrospinal fluid (CSF)

Hands-On Time: approximately 2 minutes

Sample Volume: 0.2 mL

Performance: 94.2% sensitivity, 99.8% specificity⁵

Storage Conditions: all kit components stored at room temperature (15-25 °C)

Part Number

BIOFIRE ME Panel Reagent Kit (30 Pouches): RFIT-ASY-0118

BIOFIRE ME Panel Reagent Kit (6 Pouches): RFIT-ASY-0119

BIOFIRE® FILMARRAY® GASTROINTESTINAL (GI) PANEL

22 Targets. ~1 Hour.

BACTERIA

Campylobacter (*C. jejuni*/*C. coli*/*C. upsaliensis*)
Clostridioides (*Clostridium*) *difficile* (toxin A/B)
Plesiomonas shigelloides
Salmonella
Vibrio (*V. parahaemolyticus*/*V. vulnificus*/*V. cholerae*)
 Vibrio cholerae
Yersinia enterocolitica
Diarrheagenic *Escherichia coli*/*Shigella*
 Enteroaggregative *E. coli* (EAEC)
 Enteropathogenic *E. coli* (EPEC)
 Enterotoxigenic *E. coli* (ETEC) *lt/st*
 Shiga-like toxin-producing *E. coli* (STEC) *stx1/stx2*
 E. coli O157
 Shigella/Enteroinvasive *E. coli* (EIEC)

VIRUSES

Adenovirus F40/41
Astrovirus
Norovirus GI/GII
Rotavirus A
Sapovirus (I, II, IV, and V)

PARASITES

Cryptosporidium
Cyclospora cayetanensis
Entamoeba histolytica
Giardia lamblia

BIOFIRE GI Panel Specifications

Sample Type: stool sample in Cary Blair

Hands-On Time: approximately 2 minutes

Sample Volume: 0.2 mL

Performance: 98.5% sensitivity, 99.2% specificity⁴

Storage Conditions: all kit components stored at room temperature (15-25 °C)


Part Number

BIOFIRE GI Panel Reagent Kit (30 Pouches): RFIT-ASY-0116

BIOFIRE GI Panel Reagent Kit (6 Pouches): RFIT-ASY-0104

**PANEL BIOFIRE®
JOINT INFECTION
CERTIFIE FDA, CE2797**

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
TEST PCR NICHEE MULTIPLEX RAPIDE

BIOFIRE® Joint Infection Panel testing

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
>1 mL
synovial fluid


Analysis
CBC w/ Diff
Crystals
Biomarkers

Culture
Aerobic
Anaerobic
+/- Blood culture
bottle inoculation

Gram stain

Maldi-TOF

BioFire® Joint Infection (JI) Panel		BIO F I R E
www.BioFireDx.com		
Run Summary		
Sample ID:	123456789	Run Date: 01 Jan 2020 12:00 AM
Organisms Detected:	Staphylococcus aureus	Controls: Passed
Applicable Antimicrobial Resistance Genes Detected:	mecA/C and MREJ (MRSA)	
<p>Note: Antimicrobial resistance can occur via multiple mechanisms. A Not Detected result for a genetic marker of antimicrobial resistance does not indicate susceptibility to associated antimicrobial drugs or drug classes. A Detected result for a genetic marker of antimicrobial resistance cannot be definitively linked to the microorganism(s) detected. Culture is required to obtain isolates for antimicrobial susceptibility testing and BioFire Joint Infection Panel results should be used in conjunction with culture results for the determination of susceptibility or resistance.</p>		
Result Summary		
Antimicrobial Resistance Genes		
☒	N/A	CTX-M
☒	N/A	IMP
☒	N/A	KPC
☑	Detected	mecA/C and MREJ (MRSA)
☒	N/A	NDM
☒	N/A	OXA-48-like
☒	N/A	vanA/B
☒	N/A	VIM
Gram Positive Bacteria		
Not Detected	<i>Anaerococcus prevotii/vaginalis</i>	
Not Detected	<i>Clostridium perfringens</i>	
Not Detected	<i>Cutibacterium avidum/granulosum</i>	
Not Detected	<i>Enterococcus faecalis</i>	
Not Detected	<i>Enterococcus faecium</i>	
Not Detected	<i>Finnegoldia magna</i>	
Not Detected	<i>Parvimonas mirra</i>	
Not Detected	<i>Peptoniphilus</i>	
Not Detected	<i>Peptostreptococcus anaerobius</i>	
Detected	<i>Staphylococcus aureus</i>	
Not Detected	<i>Staphylococcus lugdunensis</i>	
Not Detected	<i>Streptococcus</i> spp.	
Not Detected	<i>Streptococcus agalactiae</i>	
Not Detected	<i>Streptococcus pneumoniae</i>	
Not Detected	<i>Streptococcus pyogenes</i>	
Gram Negative Bacteria		
Not Detected	<i>Bacteroides fragilis</i>	
Not Detected	<i>Citrobacter</i>	
Not Detected	<i>Enterobacter cloacae</i> complex	
Not Detected	<i>Escherichia coli</i>	
Not Detected	<i>Haemophilus influenzae</i>	
Not Detected	<i>Kingella kingae</i>	
Not Detected	<i>Klebsiella aerogenes</i>	
Not Detected	<i>Klebsiella pneumoniae</i> group	
Not Detected	<i>Morganella morganii</i>	
Not Detected	<i>Neisseria gonorrhoeae</i>	
Not Detected	<i>Proteus</i> spp.	
Not Detected	<i>Pseudomonas aeruginosa</i>	
Not Detected	<i>Salmonella</i> spp.	
Not Detected	<i>Serratia marcescens</i>	
Yeast		
Not Detected	<i>Candida</i>	
Not Detected	<i>Candida albicans</i>	
Run Details		
Pouch:	JI Panel v1.1	Protocol: SF v3.1
Run Status:	Completed	Operator: Anonymous
Serial No.:	01234567	Instrument: FA0000
Lot No.:	012345	

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1. Overall performance based on prospective clinical study for the BIOFIRE® Joint Infection Panel, data on file, BIOFIRE Diagnostics

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


Le panel BIOFIRE® Joint Infection (JI)

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
Co-detections by organism


1. Tande JT, Patel R, Clin Microbiol Rev. 2014 Apr; 27(2): 302–345.

2. Graue C, et al. Evaluation of the BioFire Bone and Joint Infection (BJI) Panel for the Detection of Microorganisms and Antimicrobial Resistance Genes in Synovial Fluid Specimens. IDWeek 2020; Virtual.

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

5-6 Weeks Time to Result!

- GP starts abx

- Orthopedist performs additional clinical exam
- Orders radiographs and examines imaging done at outside institution

- Draws blood and orders ESR, CRP, WBC
- Aspirates joint and orders neutrophil count and differential, gram stain, and culture


- Halts abx, orders reaspiration in 3 weeks

- Synovial fluid WBC + diff consistent with PJI
- Gram stain negative
- Culture positive for *Staphylococcus aureus* after 3 days
- Surgical treatment planned

Results could potentially be available 4 weeks earlier!

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
- Physical exam
- CRP and ESR blood tests ordered
- X-ray performed

- Grossly purulent, bloody joint fluid collected
- WBC count and differential performed

- Gram stain of fluid negative for microorganisms
- Fluid inoculated on solid media and blood culture bottle

- Gram stain from bottle shows short, plump Gram-neg rods
- Positive bottle broth subcultured

- Common colonizer of young children
- Can hematogenously spread
- A common cause of septic arthritis and osteomyelitis in children 6-36 months
- Facultative anaerobe, very fastidious

Results could potentially be available 4-5 days earlier!

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MERCI INFINIMENT!

Scan the QR code
for more information.

Learn more about the BIOFIRE range of commercially-available panels for syndromic infectious disease diagnostics.

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PIONEERING DIAGNOSTICS